





Product Advantages

- Spot diameter and interval are adjustable
- The best USA RF laser tube


Continuous and super pulse

- Easy and fast treatment
- For different area of the body










Imported crystal bonding laser rod, best quality, shot more than 1.5 million times Counter on the handle, more easy to get the total shot time

LED red indicator, more easy to aim on the treatment area. 5 treatment tips for multi fuctions: tattoo removal, skin


rejuvenation, pigment removal, carbon peeling....

Treatment Effect









Screen Image

ND YAG LASER SYSTEM

Ω Skin Rejuvenation Nevus of Ota 8.4 inch true color touch screen. Super cool touch Experience. Special language can be added according to your actual requirements. Your logo also can be added

Weight

Package Size

Details of machine

Parameters Table



on diode laser so	reen.
Product name	Q switch NDYAG Laser Machine
Model	TELEVA ND YAG
Probe	Adjustable 1064/532nm probe(1-9mm spot size) . Fixed 1064/532nm probe . 1320nm probe for carbon peeling treatment
Wavelength	1064&532&1320nm
Energy	1-1500MJ
Width of Pulse	0-10 ns
Frequency	1-10 Hz
Power Supply	500W Real Work power

N.W: 22KG; G.W: 30KG

61*48*48 CM